

# **Cool**Slimming

WELCOME

### Freeze Stubborn Fat away, Say NO to surgery


The state of-the-art freezing technology removes excess fat from a specific part of the body by dissolving fat cells using cool temperature without damaging the surrounding tissue. Once cold, fat cells will decompose naturally and gently be eliminated through the body's normal metabolic processes.



ECM

ŢŗįÇęļļ

TUN







#### Sit back, Relax and Say Goodbye to Stubborn Fat

The exclusive design of Cool Slimming enables us to work with four applicators simultaneously.

- Shorter treatment time
- Greater patient comfort
- Homogenous effect
- Visible results in 30 days

#### **Application:**

- Localized fat reduction
- Stubborn fat elimination
- Target areas













Chin

Waist

Arms

Thighs


Abdomen

#### **Features:**

- Freezing temperature range: -11  $^{\circ}$  C  $\sim$  5  $^{\circ}$  C
- Enhanced cup design offering 360° freezing.
- 60% more areas treated.
- Visible results in 30 days.
- Three different cups sizes for multiple target areas.
- New patented four handle technology + chin handle
- Includes Massaging + pre heating
- Multi-patient modality



80\*40MM MINI APPLICATOR FOR DOUBLE CHIN TREATMENT





# Cool Slimming 360° - Application Parameter Setting

80/130/160/220	Thickness   <20 mm   20-40 mm   40-50 mm   50-60 mm   >60 mm	Dose 0°C -5°C -8°C -10°C -11°C	<b>Time</b> 60 min 45 min 45 min 45 min 45 min
80/130/160/220 Solution Note: Please operate setting is under -5 de	-	Dose -4°C -7°C -8°C -10°C	Time 60 min 45 min 45 min 40 min

## **Results:**



11



Most data presented in this brochure are based on published studies carried out in different parts of the world. Health-related information and opinions change frequently and therefore information contained in this brochure may be outdated, incomplete or incorrect. The information, statements and images provided in this brochure is for informational purposes only and should not substitute advice provided by the treating physician.

Manufactured by: PZ LASER SLIM TECHNOLOGY CO.LIMITED #Room 201, Floor 2nd, Buliding 2, No.52 Hongsong Road, High-Tech Industrial Development Zone, Zhengzhou, R.O.C CE, TUV, FDA and ECM are registered with PZ Laser Slim Technology Co. Limited