

# S-CURVE LIPOLASER

An expert solution to shape  
and contour your body!!



**MELT YOUR FAT DOWN !!**

S-Curve is the latest in non-invasive fat reduction, offering you a way to achieve selective lipolysis at specified areas and body contouring with

**NO PAIN,  
NO NEEDLES and  
NO DOWNTIME.**


## An expert solution to shape and contour your body

The S-Curve Lipolaser emits low levels of visible red laser energy, which creates a chemical signal in the fat cells, breaking down the stored triglycerides into fatty acids and glycerol thereby releasing them through channels in the cell membranes. The fatty acids and glycerol are metabolized and then transported around the body to create energy.


This is a naturally occurring process, thus S-Curve Lipolaser is not creating any unnatural reaction in the body nor does it affect or damage any surrounding structures such as skin, blood vessels and peripheral nerves. A period of exercise post treatment will ensure the complete metabolism and elimination of the free fatty acids from the body .

### Mechanism of Action


- Operates on Dual wavelength technology
- Bio stimulation effect on fat tissues.
- **650 nm** penetrates subcutaneous layer and give fat dissolving.
- **940 nm** penetrates deeper and accelerate lymphatic drainage, blood circulation and metabolism.




*Fat cell rest status*




*Laser irradiate fat cell*




*Penetrate fat cell, outflow the intracellular material*



*Water, glycerol, free fatty acid been transferred to mesenchyme space*



*More water (blue), free fatty acid (green) and glycerol leakage*



*Fat cell shrink significantly*

### Innovative Paddles

- Intuitively designed laser paddles with inbuilt Mitsubishi Diodes.
- 130\*60mm, 1.5 times larger than similar products.
- Customized diodes and paddles for multiple target areas.



## Features

### ▶ Clinically Proven

Ultrasound images shows up to 30% reduction in the fat layer depth after just one treatment.


### ▶ Affordable Treatment

Compared to surgical liposuction and other Ultrasound or laser techniques, the S Curve Lipo Laser is much more affordable while offering almost similar results.


### ▶ Course Treatment

A course of 8 treatments is recommended over 4 weeks with 2 treatments per week.


### ▶ Visible results in just few sessions.




## Applications




Chin




Waist



Arms



Love Handles



Thighs

## Post Treatment

Follow up with cardio and exercise for faster metabolism is mandatory.

## Results

Typically 1-2 cm loss in abdomen or target areas with visible results in a few sessions.



BEFORE



AFTER



BEFORE



AFTER



BEFORE



AFTER



BEFORE



AFTER

## Testimonials



*I am confident that I can offer my patients a treatment that is going to improve their cellulite. In fact, S-Curve Lipolaser is a great investment for building patient loyalty, while improving our medical prestige.*

**- From Austria Dr Ben**



*S-Curve Lipolaser is my workhorse body shaping product. It is a must have for any cosmetic or dermatologic surgeon offering surgical and non surgical body shaping.*

**- USA Dr Ana**

Brought to you exclusively by:



TriCell Biologics Pvt. Ltd  
Plot No: 4&5, Dhanarajpuram, Kolathur,  
Chennai - 600 099, Tamilnadu, India.  
Tel. +91 95000 80399  
E-mail: support@tricellbio.com  
Website: www.tricellbio.com



### Disclaimer:

Most data presented in this brochure are based on published studies carried out in different parts of the world. Health-related information and opinions change frequently and therefore information contained in this brochure may be outdated, incomplete or incorrect. The information, statements and images provided in this brochure is for informational purposes only and should not substitute advice provided by the treating physician.