

COOL SLIMMING 360°

TECHNICAL PARAMETERS

- (1) Rated input power: AC110V/60Hz (PZ-816A) or AC220V/50Hz (PZ-816B)
- (2) Maximum input power: 1500W
- (3) Heating temperature range: 37 ° C ~ 42 ° C
- (4) Refrigeration temperature range: -11 ° C ~ 5 ° C
- (5) Negative pressure: 10~80Kpa
- (6) Treatment handle indicator: red (630nm) / green (570nm) / 50mW × 4
- (7) Coolant: well filtered deionized or distilled water
- (8) Volume: 130cm × 62cm × 62cm
- (9) Total weight (net weight): 75kg
- (10) Operation ability: can work continuously for 24 hours
- (11) Fuse type: F250V L20A (PZ-816A) or F250V L10A (PZ-816B)
- (12) Working temperature: 10 ° C ~ 40 ° C
Storage temperature: 5 ° C ~ 55 ° C
- (13) Relative humidity: ≤80%
- (14) Atmospheric pressure: 86Kpa~106Kpa