

Eos-Ice Plus

The golden standard for laser hair removal



USPs of the Device

- 3-in-1 wavelength Technology



- Targets multiple areas of Hair follicle


- Twin Cooling System
 - Suitable for all skin type
 - Half the treatment time
- Extra Cooler handle provided for patient comfort



Technical Parameters	Device Specifications
Laser Wavelength	755 nm/808nm/1064 nm
Pulse Frequency Range	1~10Hz , 15 Hz, 20 Hz
Pulse Width Range	10-400 ms
Energy Range	1~ 120 J/cm ²
Cooler System	Twin Cooling System : Intregrated Sapphire Tip(-3°C) + Cooler Handle
Refrigeration system	Water+Semiconductor+Air cooling
Handle pieces	Laser Handle +Cooler Handle
Maximum Input power	2000W
Laser Life Span	20, 000,000 million shots
Rated Input power	AC 110V/60Hz AC220V/50Hz
Working modes	Expert and Intelligent
Total weight	85 kg
Volume	132 cm x 63xm x 57 cm
Coolant	Well-filtered distilled or deionized water
Operation temperature	10°C ~ 32°C
Storage Temperature	5°C ~ 55°C
Treatment Handle Indicator	Green
Operation Capability	Supports 24 Hours of Uninterrupted Work




Applications




Chin




Back



Beard



Leg



Chest



Lip



Bikini



Arm